

# CAPELLA



## *Elina M3* Foetal Monitoring System

### FEATRURES

- 12.1 inch true LCD flip screen
- Hidden handle, easy to move
- Low ultrasound probe with high sensitive multi-Wafer(1MHz)
- FHR 120- 160bpm normal area mark
- Automatic data storage of 48 hours for play back or print out.
- Built-in long life-span thermal printer, MTBF more than 20 years, print immediately or after storing with clear and exquisite curve.
- 152mm wide paper, support American and international standard paper.
- Setting up printing time, automatic stopping printing at appointed time.
- DSP tech, more advanced and convenient. With event marker.
- Code-operation, easy and convenient.
- A variety of alarming functions, include tachycardia, bradycardia, printing time-over indicator, shortage of paper, etc.
- With manual fetal movement functions.
- Standard configuration: single FHR monitoring



### OPTIONAL

- Twin FHR monitoring.
- Maternal SpO2, HR, NIBP, Respiration Rate, ECG, Temperature

## CAPELLA

Capella Electronics Limited  
Unit 1213, Floor 12, Wing On Plaza  
62 Mody Road, Tsim Sha Tsui, Hong Kong  
Tel: (852) 81229497 / 31668514 Fax: (852) 31668414  
URL: [www.capellaelectronics.com.hk](http://www.capellaelectronics.com.hk)  
Email: [info@capellaelectronics.com.hk](mailto:info@capellaelectronics.com.hk)  
[sales@capellaelectronics.com.hk](mailto:sales@capellaelectronics.com.hk)

\*Appearance and specifications may subject to change without any prior notice

CE ISO 13485  
certified

© Capella Electronics Limited